



TRAP240 rabbit pAb antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (kD):
A22831	Rabbit	1 mg/ml	239 kD
Applications	IHC,WB,ELISA		
Reactivity	Human		
Dilution	WB 1:500-2000 IHC: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.		
Storage	-20°C/1 year		
Specificity	TRAP230 Polyclonal Antibody detects endogenous levels of TRAP230 protein.		
Source / Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.		
Immunogen	The antiserum was produced against synthesized peptide derived from human MED12. AA range:611-660		
Uniprot No	Q93074		
Alternative names	MED12; ARC240; CAGH45; HOPA; KIAA0192; TNRC11; TRAP230; Mediator of RNA polymerase II transcription subunit 12; Activator-recruited cofactor 240 kDa component; ARC240; CAG repeat protein 45; Mediator complex subunit 12; OPA-containing prote		
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.		
Clonality	Polyclonal		
Isotype			
Conjugation			
Background	mediator complex subunit 12(MED12) Homo sapiens The initiation of transcription is controlled in part by a large protein assembly known as the preinitiation complex. A component of this preinitiation complex is a 1.2 MDa protein aggregate called M		
Other	Gene_name: MED12 ; Protein_name: Mediator of RNA polymerase II transcription subunit 12; Expression: Bone marrow,Brain,Cervix carcinoma,Epithelium,		
Product Images	<input type="checkbox"/>		

Application Key:



W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

Species Cross-Reactivity Key:

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.