



CstF-50 rabbit pAb antibody

| Catalog No : | Source: | Concentration : | Mol.Wt. (kD): |
|--------------|---------|-----------------|---------------|
| A13136 | Rabbit | 1 mg/ml | 48 kD |

| | |
|------------------------------|---|
| Applications | WB,IHC,ELISA |
| Reactivity | Human,Mouse,Rat |
| Dilution | WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications. |
| Storage | -20°C/1 year |
| Specificity | CstF-50 Polyclonal Antibody detects endogenous levels of CstF-50 protein. |
| Source / Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen | The antiserum was produced against synthesized peptide derived from human CSTF1. AA range:1-50 |
| Uniprot No | Q05048 |
| Alternative names | CSTF1; Cleavage stimulation factor subunit 1; CF-1 50 kDa subunit; Cleavage stimulation factor 50 kDa subunit; CSTF 50 kDa subunit; CstF-50 |
| Form | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Clonality | Polyclonal |
| Isotype | |
| Conjugation | |
| Background | cleavage stimulation factor subunit 1(CSTF1) Homo sapiens This gene encodes one of three subunits which combine to form cleavage stimulation factor (CSTF). CSTF is involved in the polyadenylation and 3' end cleavage of pre-mRNAs. Similar to mamm |
| Other | Gene_name: CSTF1 ; Protein_name: Cleavage stimulation factor subunit 1; Expression: Muscle,Testis, |

Product Images

Application Key:

W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

**Species Cross-Reactivity Key:**

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.