



Ku-80 rabbit pAb antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (kD):
A16859	Rabbit	1 mg/ml	80 kD
Applications	WB,IHC,IF,ELISA		
Reactivity	Human,Monkey		
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. IF: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.		
Storage	-20°C/1 year		
Specificity	Ku-80 Polyclonal Antibody detects endogenous levels of Ku-80 protein.		
Source / Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.		
Immunogen	The antiserum was produced against synthesized peptide derived from human Ku70/80. AA range:683-732		
Uniprot No	P13010		
Alternative names	XRCC5; G22P2; X-ray repair cross-complementing protein 5; 86 kDa subunit of Ku antigen; ATP-dependent DNA helicase 2 subunit 2; ATP-dependent DNA helicase II 80 kDa subunit; CTC box-binding factor 85 kDa subunit; CTC85; CTCBF; DNA repair pr		
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.		
Clonality	Polyclonal		
Isotype			
Conjugation			
Background	X-ray repair cross complementing 5(XRCC5) Homo sapiens The protein encoded by this gene is the 80-kilodalton subunit of the Ku heterodimer protein which is also known as ATP-dependant DNA helicase II or DNA repair protein XRCC5. Ku is the DNA-bin		
Other	Gene_name: XRCC5 ; Protein_name: X-ray repair cross-complementing protein 5; Expression: Cervix carcinoma,Coronary artery,Heart,Neuroblastoma,Osteoblast,Thy		

Product Images

Application Key:



W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

Species Cross-Reactivity Key:

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.