



Neurensin-1 rabbit pAb antibody

| Catalog No : | Source: | Concentration : | Mol.Wt. (kD): |
|--------------|---------|-----------------|---------------|
| A18257 | Rabbit | 1 mg/ml | 25 kD |

| | |
|-----------------------|---|
| Applications | WB,IHC,ELISA |
| Reactivity | Human,Mouse,Rat |
| Dilution | WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications. |
| Storage | -20°C/1 year |
| Specificity | Neurensin-1 Polyclonal Antibody detects endogenous levels of Neurensin-1 protein. |
| Source / Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen | The antiserum was produced against synthesized peptide derived from human NRSN1. AA range:121-170 |
| Uniprot No | Q8IZ57 |
| Alternative names | NRSN1; VMP; Neurensin-1; Neuro-p24; Vesicular membrane protein of 24 kDa; Vesicular membrane protein p24 |
| Form | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Clonality | Polyclonal |
| Isotype | |
| Conjugation | |
| Background | function:May play an important role in neural organelle transport, and in transduction of nerve signals or in nerve growth. May play a role in neurite extension.,similarity:Belongs to the VMP family.,subcellular location:Localizes mainly to neurites.,tiss |
| Other | Gene_name: NRSN1 ; Protein_name: Neurensin-1; Expression: Brain,Fetal brain,Teratocarcinoma, |

Product Images

Application Key:

W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

**Species Cross-Reactivity Key:**

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.