



Olfactory receptor 5L1/2 rabbit pAb antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (kD):
A18959	Rabbit	1 mg/ml	35 kD
Applications	WB,IF,ELISA		
Reactivity	Human		
Dilution	WB: 1/500 - 1/2000. IF: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.		
Storage	-20°C/1 year		
Specificity	Olfactory receptor 5K3 Polyclonal Antibody detects endogenous levels of Olfactory receptor 5K3 protein.		
Source / Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.		
Immunogen	The antiserum was produced against synthesized peptide derived from human OR5K3. AA range:232-281		
Uniprot No	A6NET4		
Alternative names	OR5K3; Olfactory receptor 5K3		
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.		
Clonality	Polyclonal		
Isotype			
Conjugation			
Background	olfactory receptor family 5 subfamily K member 3(OR5K3) Homo sapiens Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members o		
Other	Gene_name: OR5K3 ; Protein_name: Olfactory receptor 5K3; Expression:		
Product Images			

Application Key:

W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

**Species Cross-Reactivity Key:**

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.