



p140Cap rabbit pAb antibody

| Catalog No : | Source: | Concentration : | Mol.Wt. (kD): |
|-----------------------|--|-----------------|---------------|
| A19232 | Rabbit | 1 mg/ml | 112 kD |
| Applications | WB,IHC,IF,ELISA | | |
| Reactivity | Human | | |
| Dilution | WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. IF: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications. | | |
| Storage | -20°C/1 year | | |
| Specificity | p14 Polyclonal Antibody detects endogenous levels of p14 protein. | | |
| Source / Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. | | |
| Immunogen | The antiserum was produced against synthesized peptide derived from human p14 ARF. AA range:71-120 | | |
| Uniprot No | Q8N726 | | |
| Alternative names | CDKN2A; CDKN2; MLM; Cyclin-dependent kinase inhibitor 2A; isoform 4; p14ARF; p19ARF | | |
| Form | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. | | |
| Clonality | Polyclonal | | |
| Isotype | | | |
| Conjugation | | | |
| Background | CDKN2A generates several transcript variants which differ in their first exons. At least three alternatively spliced variants encoding distinct proteins have been reported, two of which encode structurally related isoforms known to function as inhibitors | | |
| Other | Gene_name: CDKN2A ; Protein_name: Cyclin-dependent kinase inhibitor 2A isoform 4; Expression: | | |
| Product Images | | | |

Application Key:

W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

**Species Cross-Reactivity Key:**

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.