



TRAP230 rabbit pAb antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (kD):
A22830	Rabbit	1 mg/ml	247 kD
Applications	WB,IF,ELISA		
Reactivity	Human,Mouse		
Dilution	WB: 1/500 - 1/2000. IF: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.		
Storage	-20°C/1 year		
Specificity	TRAP220 Polyclonal Antibody detects endogenous levels of TRAP220 protein.		
Source / Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.		
Immunogen	The antiserum was produced against synthesized peptide derived from human PPAR-BP. AA range:1423-1472		
Uniprot No	Q15648		
Alternative names	MED1; ARC205; CRSP1; CRSP200; DRIP205; DRIP230; PBP; PPARBP; PPARGBP; RB18A; TRAP220; TRIP2; Mediator of RNA polymerase II transcription subunit 1; Activator-recruited cofactor 205 kDa component; ARC205; Mediator complex subunit 1; Peroxiso		
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.		
Clonality	Polyclonal		
Isotype			
Conjugation			
Background	mediator complex subunit 1(MED1) Homo sapiens The activation of gene transcription is a multistep process that is triggered by factors that recognize transcriptional enhancer sites in DNA. These factors work with co-activators to direct transcri		
Other	Gene_name: MED1 ; Protein_name: Mediator of RNA polymerase II transcription subunit 1; Expression: Cervix carcinoma,Colon,Epithelium,Heart,Testis,		
Product Images			

Application Key:



W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

Species Cross-Reactivity Key:

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.