



## C9 rabbit pAb antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (kD):
A11428	Rabbit	1 mg/ml	70 kD
Applications	WB,ELISA		
Reactivity	Human		
Dilution	WB: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.		
Storage	-20°C/1 year		
Specificity	C9 Polyclonal Antibody detects endogenous levels of C9 protein.		
Source / Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.		
Immunogen	The antiserum was produced against synthesized peptide derived from human C9. AA range:181-230		
Uniprot No	P02748		
Alternative names	C9; Complement component C9		
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.		
Clonality	Polyclonal		
Isotype			
Conjugation			
Background	complement C9(C9) Homo sapiens This gene encodes the final component of the complement system. It participates in the formation of the Membrane Attack Complex (MAC). The MAC assembles on bacterial membranes to form a pore, permitting disruption o		
Other	Gene_name: C9 ; Protein_name: Complement component C9; Expression: Blood,Liver,Plasma,Platelet,		

### Product Images

#### Application Key:

W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation  
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

#### Species Cross-Reactivity Key:



H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster  
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All  
Species Expected

**Trademarks**

*All product names and trademarks are the property of their respective owners.*

**Regulatory Disclaimer**

*For life science research only. Not for use in diagnostic procedures.*

---

**Contact and Support:**

*To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).*

*To call, write, fax, or email us, please visit [www.aabsci.com](http://www.aabsci.com), contact information will be displayed.*