



## CD24 rabbit pAb antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (kD):
A11959	Rabbit	1 mg/ml	kD
Applications	IHC,ELISA		
Reactivity	Human		
Dilution	IHC 1:50-200, ELISA 1:10000-20000		
Storage	-20°C/1 year		
Specificity	The antibody detects endogenous CD24		
Source / Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.		
Immunogen	Synthetic peptide from human protein at AA range: 30-60		
Uniprot No	P25063		
Alternative names	Signal transducer CD24 (Small cell lung carcinoma cluster 4 antigen) (CD antigen CD24)		
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.		
Clonality	Polyclonal		
Isotype			
Conjugation			
Background	CD24 molecule(CD24) Homo sapiens This gene encodes a sialoglycoprotein that is expressed on mature granulocytes and B cells and modulates growth and differentiation signals to these cells. The precursor protein is cleaved to a short 32 amino acid m		
Other	Gene_name: CD24 CD24A ; Protein_name: Signal transducer CD24 (Small cell lung carcinoma cluster 4 antigen) (CD antigen CD24); Expression: Ovary,		
Product Images			

### Application Key:

W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation  
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

### Species Cross-Reactivity Key:



H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster  
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All  
Species Expected

**Trademarks**

*All product names and trademarks are the property of their respective owners.*

**Regulatory Disclaimer**

*For life science research only. Not for use in diagnostic procedures.*

---

**Contact and Support:**

*To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).*

*To call, write, fax, or email us, please visit [www.aabsci.com](http://www.aabsci.com), contact information will be displayed.*